

ONE-HAND OPERATED TURBO JOYSTICK FOR NINTENDO

--- DISCONTINUED. (Verified 05/2008) RETAINED IN DATABASE FOR REFERENCE. ---

The One-Hand Operated Turbo Joystick for Nintendo is a joystick designed to enable persons with the use of a hand or the mouth to operate the Nintendo and Genesis systems. The joystick is mounted on a base which can be secured to a table. A flexible stand holds the sip/puff tube in an accessible position. A turbo switch gives rapid fire capability for A and B buttons by means of pneumatic controls, and a lever for the Start and Select functions. The units plug directly into the Nintendo game system without modification and have a standard Nintendo plug so that a friend can join in and help with difficult game situations.

Notes: Adapter cables for both joystick models are available for Genesis (\$29), Sega (\$10), or Atari (\$10).

Price: Contact manufacturer.

This product record was updated on January 18, 2001.

This product is available from:

Manufacturer:

KYE (Formerly KY Enterprises)

12 Mathew Bird Circle

Bozeman, Montana 59715

United States

Telephone: 406-586-2376.

Fax: 646/304-6904.

Web: <http://www.quadcontrol.com>.

Email: KenYankelevitz@msn.com.

TWO-HAND-OPERATED JOYSTICK FOR NINTENDO

--- DISCONTINUED. (Verified 05/2008) RETAINED IN DATABASE FOR REFERENCE.

The Two-Hand-Operated Joystick for Nintendo is designed to allow people with limited use of both hands to operate the Nintendo games. The device is a modified NES Advantage joystick with turbo rapid fire and remote A and B buttons.

Notes: Shipping & handling and California sales tax not included in price listed.

Price: Contact manufacturer.

This product record was updated on May 21, 2008.

This product is available from:

Manufacturer:

KYE (Formerly KY Enterprises)

12 Mathew Bird Circle

Bozeman, Montana 59715

United States

Telephone: 406-586-2376.

Fax: 646/304-6904.

Web: <http://www.quadcontrol.com>.

Email: KenYankelevitz@msn.com.

[Link to more products from KYE \(Formerly KY Enterprises\)](#)

AbleData does not produce, distribute or sell any of the products listed on this website, but we provide you with information on how to contact manufacturers or distributors of these products. If you are interested in purchasing a product, you can find companies who sell it below.

Ezra:

DISCONTINUED. (Verified 05/2008) RETAINED IN DATABASE FOR REFERENCE.

EZRA is an environmental control system that can be operated by a single switch, a puff/sip (pneumatic) switch, or a button switch. The EZRA base unit (\$950) can control television and radio (radio is built-in), turn lights and appliances on and off, control a VCR, and unlock an electric door. The unit can also control an electric bed, control an intercom system, make a phone call, control a compact disk player, and access a 50-word phrase list for non-verbal individuals with additional units. Another option is to use a Universal Remote (\$100), which does not require having the TV, VCR, and cable remotes modified (these will be modified at no additional cost). The system is menu driven, with the list of options displayed on a TV screen with an arrow that constantly moves down the list at the speed selected by the user. POWER: 120 volt AC. DIMENSIONS: 19 x 4 x 14 inches. WARRANTY: 90-day warranty; 30-day moneyback guarantee (less shipping costs).

Price 501 To 1000 Dollars

Price 1001 To 5000 Dollars

Fingertip Controller

DISCONTINUED. (Verified 06/93) RETAINED IN DATABASE FOR REFERENCE.

FingerTip Controller includes five soft-touch buttons which simulate the action of a standard joystick to control screen movement in four directions, with the fifth button to fire or initiate action. A left hand model of the controller is available. The FingerTip Controller is connected to the computers with a 9-pin DIN jack. COMPATIBILITY: For use on Apple II, Atari (some or all models), Commodore 64, and NEC (some or all models) computers.

Price Check

Price: 15.00 as of: 10/02/1997

Mouth Operated Controller For Quads

AbleData does not produce, distribute or sell any of the products listed on this website, but we provide you with information on how to contact manufacturers or distributors of these products. If you are interested in purchasing a product, you can find companies who sell it below.

DISCONTINUED. (Verified 06/93) RETAINED IN DATABASE FOR REFERENCE.

Mouth Operated Controller for Quads is an alternative joystick that requires only 1 ounce of pressure for activation. The user moves a lever with his tongue in order to execute joystick movements and sips and puffs on a tube in order to operate the fire, reset, and game select buttons. The system functions as a directional control that keeps rate constant and can change direction. The system is fully compatible with the Atari computers; some modifications must be made for some of the other computers. These can be done by KY Enterprises for the cost of postage, or KY Enterprises will send instructions. Contact manufacturer for further information. COMPATIBILITY: For use on Apple II, Atari (some or all models), and Commodore (some or all models) computers.

Price Check

Price: 50.00 as of: 10/02/1997

One-Hand Operated Turbo Joystick For Nintendo

DISCONTINUED. (Verified 05/2008) RETAINED IN DATABASE FOR REFERENCE.

The One-Hand Operated Turbo Joystick for Nintendo is a joystick designed to enable persons with the use of a hand or the mouth to operate the Nintendo and Genesis systems. The joystick is mounted on a base which can be secured to a table. A flexible stand holds the sip/puff tube in an accessible position. A turbo switch gives rapid fire capability for A and B buttons by means of pneumatic controls, and a lever for the Start and Select functions. The units plug directly into the Nintendo game system without modification and have a standard Nintendo plug so that a friend can join in and help with difficult game situations.

Portable Self Contained Video Game Unit

DISCONTINUED. (Verified 06/93) RETAINED IN DATABASE FOR REFERENCE. Video game unit integrates Atari video game system, TV and controls for joystick, paddle, game select and reset, and game difficulty levels.

400 up Dollars (Price check 8/4/2002)

Turbo Sip/puff Joystick For Nintendo

DISCONTINUED. (Verified 05/2008) RETAINED IN DATABASE FOR REFERENCE.

The Turbo Sip/Puff Joystick for Nintendo is a modified game controller that allows the joystick to be mouth operated by people with disabilities. To use, the user's mouth is placed over a lollypop shaped plastic disk. Head movement is required to operate the joystick and to push the start/select lever switch. The Turbo and reset switches are included. Turbo gives the user rapid fire mode for the A and B switches. To make the reset switch work, two wires to the circuit board of the Nintendo game must be soldered. The joystick is mounted on a stand which can be secured to a table or wheelchair tray. Instructions are included. DIMENSIONS: 4 x 4 inches. WEIGHT: 1 pound.

Tongue Controller

DISCONTINUED. (Verified 06/93) RETAINED IN DATABASE FOR REFERENCE.

Tongue controller is a modified telephone touch pad which includes four directional switches to control screen movement, and action, reset, and select switches. The Tongue Controller connects to the computer with a 9-pin DIN jack, and requires an activation pressure of 0.5 oz. to trigger the switch. COMPATIBILITY: For use on Apple II, Atari (some or all models), Commodore 64, and NEC (some or all models) computers.

Price Check

Price: 40.00 as of: 10/02/1997

Two-Hand-Operated Joystick For Nintendo

DISCONTINUED. (Verified 05/2008) RETAINED IN DATABASE FOR REFERENCE.

The Two-Hand-Operated Joystick for Nintendo is designed to allow people with limited use of both hands to operate the Nintendo games. The device is a modified NES Advantage joystick with turbo rapid fire and remote A and B buttons.

Video Game Controllers.

DISCONTINUED. (Verified 06/93) RETAINED IN DATABASE FOR REFERENCE.

Mouth operated video game controllers for use by quadriplegics. Interface with the Atari video game system. Operates game select, reset, paddle, joystick and game difficulty switches. Several models available for mouth operation: soft touch buttons, puff switches, joysticks. Mounting stand brings controls to correct height. Hand/wrist controls available. Multi-game board allows remote selection of cartridges, pushbutton switch changes the game

Price 150 Dollars

[See more on Ken Yankelevitz at: <https://www.oneswitch.org.uk/art.php?id=394>]



AbleData

Your source for assistive technology information

Alternate Layout | My AbleData | Contact | Site Map | A-Z Index | Help |

Home | Products | Resources | Library | About AbleData

"AbleData provides objective information about assistive technology products and rehabilitation equipment"



Highlights

[100 most recently added products.](#)

The [10th International Conference on Low Vision](#)- February 20, 2011 to February 24, 2011 in Kuala Lumpur, Malaysia.

The Shepherd Center in Atlanta is conducting research on the use and usability of wireless technology by people with disabilities. They need volunteers to take an online survey before December 31, 2010. The survey is brief and you will be entered in a drawing to win a gift certificate. [More Details.](#)

Follow us on Twitter to find out about new products as soon as we add them.



[Our most recent product tweets.](#)

About AbleData

- We do not sell products, but can help you locate companies that do
- Domestic and international sources

[More about AbleData](#)

Features

- [Almost 40,000 product listings in 20 categories.](#)
- [Resource listings including information centers, conferences, companies by state, and help wanted.](#)
- [Library of publications, literature, and news.](#)

MyAbleData

Customize your experience by saving products, searches, and preferences. Add product reviews, classifieds, and more.

[Sign up or login](#)

Main Menu

Hyper-ABLEDATA allows you to search for product or company information in different ways.

Choose the type of search you would like to do ---->

Search by TYPE of product

Search by PRODUCT NAME

Search by COMPANY NAME

Allows you to find products if you know part of Product name



Startup Card


Write Us


QUIT